

	example REF	Coupling	Drill head	- Size
		28440	S	- 105
Coupling	AO	25440		
	Harris	26440		
	Hudson	27440		
	Stryker/Zimmer	28440		
	Aesculap	29400		
Drill head	sharp head		S	
	round head		R	
Size Ø	6.0 mm			- 060

	22.0 mm			- 220

Medullary Reamer, flexible, Ø 10.5 mm
sharp head total length 480 mm
Stryker/Zimmer-coupling cannulated Ø 3.5 mm

size Ø 6.0 - 22.0 mm available

Coupling	AO	25455N	flexible shank, nitinol
Drill head	sharp head	91080S	Ø 8.0 mm
		91085S	Ø 8.5 mm
	
		91170S	Ø 17.0 mm
		91175S	Ø 17.5 mm
Guide wire		88007-600-030	Ø 3.0 mm, length 600 mm
		88007-900-030	Ø 3.0 mm, length 900 mm

	example REF	Coupling	- Size	Front
		20480K	- 125	- 200
Coupling	AO	20240K		
Size	Ø 8.0 mm		- 080	
	
	Ø 20.0 mm		- 200	
Front	round, 200 mm			- 200
	cutting, 100 mm			- 100

Medullary Reamer, stiff, Ø 12.5 mm, AO-coupling, total length 480 mm, twist length 200 mm, round cutting, cannulated Ø 3.5 mm

Stiff reamer, length 480 mm

size Ø 8.0 - 20.0 mm available

Die **risa GmbH** ist Ihr kompetenter Partner, wenn es um chirurgische Medizinprodukte geht.

Neben unseren Bohrern und Sägblättern fertigen wir eine Vielzahl an hochpräzisen Dreh- und Frästeilen aus verschiedenen Metallen und Kunststoffen für die Industrie.

risa is your qualified partner for surgical cutting tools in medical applications.

We are experts in manufacturing drills, reamers and saw blades as well as high precision milled and turned industrial components.



MARKRAUMBOHRER

MEDULLARY REAMER

risa fertigt Markraumborher mit starrer und flexibler Welle und Markraumborher mit Wechselkopfsystem. Im Sortiment sind Borher mit Anschlüssen für alle handelsüblichen Antriebssysteme. Bohrköpfe sind in den Größen von Ø 6.0 mm bis Ø 22.0 mm in Stufen von 0,5 mm Schritten erhältlich. Die Anfertigung von Markraumborhern in verschiedenen Längen ist möglich. Darüber hinaus bieten wir Sonderanfertigungen nach Ihren Wünschen an.

risa manufactures reamers with stiff and flexible shank as well as reamers with replaceable head systems. We offer reamers with couplings for all well-known power systems. Drill heads are available in sizes from Ø 6.0 mm bis Ø 22.0 mm with increments of 0.5 mm. Fabrication of reamers in different lengths is possible and we can offer custom-made products according to your drawings.



risa GmbH

Industriestraße 7
78234 Engen | Germany

T +49 (0)7733 9423-0
F +49 (0)7733 9423-44

info@risa-germany.eu
www.risa-germany.eu

risa GmbH

Industriestraße 7
78234 Engen | Germany

T +49 (0)7733 9423-0
F +49 (0)7733 9423-20

info@risa-germany.eu
www.risa-germany.eu

FLEXIBEL | FLEXIBLE

KANNÜLIERT | CANNULATED

ANSCHLUSS | COUPLING



Länge | length 480 mm
kannüliert | cannulated Ø 3.5 mm



AO	REF 25440...	Ø 6.0 - 22.0 mm
Harris	REF 26440...	Ø 6.0 - 22.0 mm
Hudson	REF 27440...	Ø 6.0 - 22.0 mm
Stryker/Zimmer	REF 28440...	Ø 6.0 - 22.0 mm
Aesculap	REF 29400...	Ø 6.0 - 22.0 mm

FLEXIBEL MIT WECHSELKOPFSYSTEM | FLEXIBLE WITH CHANGEABLE HEADS



Länge | length 480 mm
kannüliert | cannulated Ø 3.2 mm

Führungsdraht | guide wire Ø 3.0 mm

Länge | length 600 mm | 900 mm

AO	REF 25455N	
	REF 88007-600-030	Länge length 600 mm
	REF 88007-900-030	Länge length 900 mm

STARR | STIFF



AO	REF 20480K-...-200	Ø 8.0 - 20.0 mm
	REF 20480K-...-100	Ø 8.0 - 20.0 mm